

TC 3000

At last the solution

for all Telangiecastasias's immediate disappearance.

In response to medical requests, F&P provides an exclusive solution combining effectiveness and innocuousness.



THE THERMO-COAGULATION OF THE IIIRD MILLENIUM

TC 3000

THE THERMO-COAGULATION OF THE IIIRD MILLENIUM

INDICATIONS

TELANGIECSTASIAS

Rubis points
Red blotches
Spider veins

TECHNIQUES

High frequency wave

ACTION


Thermocoagulation

EFFECT

Immediate results
No contra indication
No side effect
No pigmentation
No loss of pigmentation

ASSETS

- Instantaneous disappearance
- Non invasive
- No scars
- 40 to 50 cm of variculas erased in only a 15 minutes session
- Easy to use
- Any kind of telangiecstasias
- Any area to treat
- Any kind of skin
- Any season
- No side effect



DIRECTIONS FOR WORKING

- Generator HF : Frequency 4MHz
- Adjustable tension from 0 to 500 V by step of 5%
- Time of emission adjustable from 0,1 to 1 second

DIMENSIONS - WEIGHT

- H : 16 cm - D : 38 cm - W : 53 cm
- 8 kg

ALIMENTATION

- Voltage : 220 to 240 Volts
- Frequency : 50 to 60 Hz
- Power : 115 VA



Parc d'Activités Fontaudin - av Descartes - 33370 Artigues-près-Bordeaux, FRANCE
Tél : + 33 (0)5 56 40 01 10 - Fax +33 (0)5 56 40 01 44 - www.formesetperformances.com
PUBLIC LTD Cie: CAPITAL 700 000 euros / IDENTIFICATION FR 96 412 076 010